

...er protein source to help get you there.
...looking to get more out of their protein supple
...ed Whey Isolate is an ultra high quality protein
...maximum support for lean muscle, strength and

UPER ADVANCED WHEY ISOLATE WORTH IT:

...er scoop with
...carbohydrates.

...Whey Isolate
...that is processed
...on to ensure an
...that contains
...and fat. These
...uses separate the
...n from non-protein
...ing a highly-purified

...further enhance mixability and
...dispersion in liquid.

- Delivers a high quality Whey Isolate protein with a high amino acid profile per gram.
- Enhanced with additional free-form amino acids to work synergistically with Whey Isolate's high concentration of aminos.
- 2 scoops contain over 9 grams of the following Branched Chain Amino Acids from protein:
Which typically provides:

- Isoleucine 2.58 g
- Leucine 4.27 g
- Valine 2.28 g


...-art manufacturing
...re used to retain the
...Protein Peptides and
...s - some other whey
...processing methods
...mycomacropetides,
...an important protein
...it.

...Whey Isolate - to

...in has been relied upon for years by athletes,
...s and anyone looking to get into their best shape. Super
...Whey Isolate provides all of the benefits Whey is known
...re by carefully removing more of the fat and lactose that
...need. By helping to quickly boost protein synthesis while
...protein breakdown, Super Advanced Whey Isolate helps
...ne most out of your workout so you don't waste any of your
...rk in and out of the gym."

...in conjunction with an intense daily exercise program and a
...balanced diet including an adequate caloric intake.

...Protein Combined with a Small Amount of Carbs Provides Greater Support for Muscle
...Strength Compared to Carbs-Only When Taken During a Six Week Training Program.



Servings Per Container / Servicios por Bolsa about 26 (149) 2 Scoops (68g) / 2Cucharas (68g) about 13

Amount Per Serving / Cantidad por porción	%Daily Value / Valor Diario		%Daily Value / Valor Diario	
Calories / Calorias	110		220	
Calories from Fat / Calorias de Grasa	5		10	
Total Fat / Total Grasa	0.5 g	1%†	1 g	2%†
Saturated Fat / Grasa Saturada	<0.5 g	1%†	0.5 g	2%†
Cholesterol / Colesterol	15 mg	5%	30 mg	9%
Total Carbohydrate / Total Carbohidratos	1 g	<1%†	2 g	1%†
Dietary Fiber / Fibra	<1 g	1%†	<1 g	2%†
Sugars / Azúcar	<1 g	††	1 g	††
Protein / Proteina	30 g	60%†	60 g	120%†
Calcium / Calcio	57 mg	6%	115 mg	11%
Phosphorus / Fosforo	46 mg	5%	93 mg	9%
Magnesium / Magnesio	12 mg	3%	25 mg	6%
Sodium / Sodio	95 mg	4%	190 mg	8%
Potassium / Potasio	70 mg	2%	140 mg	4%
Creatine / Creatina	3 g	2%	6 g	4%

†Percent Daily Values are based on a 2,000 calorie diet.
††Daily Value not established.

Ingredientes: Proteína Isolada Hidrolizada de Suero/ Hydrolyzed Whey Protein Isolate, Aminoácidos Recuperadores-Super Recovery Blend (Glycine, Taurine, L-Threonine, L-Glutamine, L-Arginine, L-Alanine, L-Lysine Hydrochloride, L-Leucine, L-Isoleucine, L-Valine), Sabores Naturales/ Natural Flavors, Cellulose Gum, Xanthan Gum, Acesulfame Potassium, Sucralose, Carrageenan, Enzyme Blend / Enzimas (Ammagen, Amylase, Protease, Cellulase, Beta-D-Galactosidase, Lipase) cookies and Cream Flavor, Creatine monohydrate.

Contains milk and wheat ingredients. / Contiene ingredientes de leche y trigo.

Typical Amino Acid Profile (milligrams per 34 g scoop**)**
Tipos de Aminoácidos en miligramos (Cuchara de 34 grs)

Essential Amino Acids / Aminoácidos Esenciales		Nonessential Amino Acids / Aminoácidos no Esenciales	
Histidine	317 mg	Alanine	547 mg
Isoleucine	1,294 mg	Arginine	2,141 mg
Leucine	2,136 mg	Aspartic Acid	872 mg
Lysine	1,867 mg	Cysteine	5,491 mg
Methionine	419 mg	Glutamic Acid	5,825 mg
Phenylalanine	588 mg	Glycine	1,393 mg
Threonine	3,308 mg	Proline	500 mg
Tryptophan***	360 mg	Serine	
Valine	1,141 mg	Tyrosine	

***L-Tryptophan is naturally occurring, not added. / L-Triptófano es de origen natural, no se añade.

EXP:12/26
LOT:07052012

Importado & Distribuido por
GLOBAL NUTRICION VITAMINAS SAS Bucaramanga, Colombia
Tel: 6998604
Carrera 28 No. 40-64 Ofc. 201A
Registro Sanitario: RSiA 161138112

...to 6-8 ounces of water o
...Serious athletes and
...consume 1-2 scoops imm

Instrucciones

Para Adultos adicionar (1) se
de agua o en alguna bebida fa
Para Atletas de alto rendimie
servicios de proteína durante e

Please note: Crystalline Taurine w
as small crystals within the powder

Super Advanced Whey Isolate i
FREE OF yeast, preservative
color.

Contraindicaciones:

Si se encuentra embarazada o enferma
algún medicamento por favor consultar c
antes de ser consumido el producto.

WARNING: If you are pregnant, nu
any medications or have any medic
consult your doctor before use.
use and consult your doctor if
reactions occur. Not intended
persons under the age of 18.

KEEP OUT OF THE REACH OF
STORE AT ROOM TEMPERATURE
CLOSED AND AVOID EXCESSIVE HEAT
YOUR PROTECTION, DO NOT USE IF
UNDER CAP IS BROKEN OR MISSING.

*These statements have not been evaluated by the Food and
Drug Administration. This product is not intended to diagnose
treat, cure or prevent any disease.



Carefully Manufactured by / Fabricado por
AMERICAN FITNESS LABORATORIES, LLC
19555 E COUNTRY ADVENTURE FL 31100

